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
Espectroscopia UV-Vis

The determination of drug tablet concentration in pharmaceutical applications for drug development using the Agilent Cary 60 UV-Vis with fiber optics

Application Note

Pharma

Author

Fyfe, DJ

Fyfe Science, West Lakes
Shore, SA 5020, Australia

Note: All of this work was conducted under controlled conditions in a GMP laboratory within a leading global pharmaceutical manufacturing company.

Summary

The Agilent Cary 60 UV-Vis spectrophotometer is the new, improved

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
Figure 1. Concentration curve for anti-malarial API standards. Calibration equation: $\text{Abs} = 0.00688 \cdot \text{Conc} - 0.01056$; correlation coefficient = 0.99976. The dissolution bath was used to facilitate the dissolution to completion, and the sample was analysed for total dissolved active ingredient after time. The concentration measurement of the unknown sample was made at the end of the dissolution reaction.

Fonte: <https://www.agilent.com/cs/library/applications/5990-7945EN.pdf>

Espectroscopia UV-Vis

Measuring the purity of low volumes of DNA at 4 °C using the Agilent Cary 60 UV-Vis spectrophotometer with fiber optics microprobe

Application Note

Pharma/Biotech

Author

*Fyfe, DJ and **Comerford, JC

*Fyfe Science, West Lakes
Shore, SA 5020, Australia

**Agilent Technologies,
679 Springvale Road,
Mulgrave 3179, Australia


Summary

The Agilent Cary 60 UV-Vis spectrophotometer is the ideal instrument for measuring small amounts of biological samples directly from their storage environment, e.g., a refrigerator. This application note demonstrates how the Cary 60 is used to measure the purity of DNA at 4 °C using the micro fiber optics accessory – resulting in significant time and cost savings without compromise in the accuracy and reproducibility of data quality.

- No more need for cuvettes
- Measure samples directly from the refrigerator at 4 °C
- Save time and money per analysis

Figure 1. The Cary 60 instrument fitted with the Fiber Optics Accessory taking a reading of 150 μ L DNA in an Eppendorf tube at 4 °C

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Simple, automated measurements of the photocatalytic properties of colorimetric species using the Agilent Cary 60 UV-Vis spectrophotometer with fiber optics

Application Note Chemicals

Author

*Fyfe, DJ and **Wang, XD

*Fyfe Science, West Lakes
Shore, SA 5020, Australia


**PFPC School of Chemistry,
University of Melbourne,
3010, Australia

Summary

The Agilent Cary 60 UV-Vis spectrophotometer is the new, improved successor to the award-winning Cary 50 UV-Vis. In this short review, this instrument platform was evaluated for its potential to measure small (40 μ L) samples of methylene blue *in situ* during exposure to high-intensity UV irradiation. Methylene blue is combined with other compounds used in a variety of applications, including use in cosmetics and sun screen products and environmental remediation in contaminated air and polluted water.

Figure 1. The Cary 60 UV-Vis fitted with the fiber optics microprobe accessory provides a simple mechanism to measure a sample *in situ* and remote to the instrument

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
Cary 60

Cary 3500

- Monitor enzymatic reactions at temperature
- Calibrate and determine sample concentration
- Perform temperature ramping experiments
- Quantify nucleotides and proteins

Cary 5000

Cary 7000 UMS


Figure 1. Optical diagram of the Cary 7000 UMS accessory

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For Research Use Only.
Not for diagnostic procedures.

Using the Agilent Cary Eclipse to measure the chemiluminescence from a ruthenium complex

Application Note

Author

Mr. Stuart Purcell and Assoc.
Professor Neil Barnett
School of Biological and
Chemical Sciences
Deakin University
Geelong, 3217, Australia

Introduction

Chemiluminescence is emission resulting from a chemical reaction, as opposed to fluorescence or phosphorescence, both of which are emission resulting from the absorption of light. Since the emission is generated without light, no source is needed when measuring the signal. In a conventional fluorescence spectrophotometer the light source must be off during the measurement phase, so as not to interfere with the signal being measured. The chemiluminescence signal is also generally weak. Achieving good signal to noise and a low limit of detection depends on how well the sample is shielded from room light. It is absolutely essential that the instrument has excellent light-sealing capabilities. The Agilent Cary Eclipse fluorescence spectrophotometer addresses both these requirements which makes it an excellent analyser for chemi- or bio- luminescent applications. With a red-sensitive detector as standard and the excellent sensitivity of the gratings in the visible region, emissions well into the red can be measured with the Eclipse, without the expense of purchasing components which allow measurement in this region.

This application note demonstrates the Cary Eclipse chemiluminescence capabilities by measuring the emission from a ruthenium complex upon reaction with codeine. Figure 1 shows a photograph of the mixing of the two non-emitting samples to produce a chemiluminescence product.

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


Figure 1: Photograph of tris(2,2' - bipyridyl)ruthenium(II) chemiluminescence via the reaction with codeine

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

Fonte: [Quantification of complex polycyclic aromatic hydrocarbons or petroleum oils in water with Cary eclipse fluorescence spectrophotometer According to astm d 5412-93 \(2000\)](#)

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Fonte: [Cytosolic expression of Green Fluorescent Protein \(GFP\) and its derivatives in the yeast *Saccharomyces cerevisiae*: Detection in vivo using the Agilent Cary Eclipse](#)

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

Figure 2. Instrumental parameters for detecting DsRed emission spectra (Emission maxima 583 nm) following excitation with 550 nm (peak excitation 558 nm)

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Examine the fluorescent properties of a common laundry detergent

The excitation and emission spectra of the powder show optical brighteners absorbing in the region between 320 and 390 nm and emitting over the 400–500 nm range. Fabrics washed in this detergent would exhibit a blue hue.

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