

AFM for biological application – Theory 2

Pieter De Beule 31 October 2022

A short introduction



Your tutors for this course & AFM for biological applications







Adelaide Miranda

Nanodevices and Nanoelectronics

- Mathematics of signal processing for scientific instrumentation
- Theory 1 13:30-15:30h 10 October 2022 UMinho Ed.1-0.18

- Theory 2 14-16h 10 October 2022 UMinho Ed.2-1.10

Topics in Biophysics and Bionanosystems

- **Bio-AFM**
 - Theory 1 15:30-17:30h
- 24 October 2022 UMinho Ed.1-1.18
- Theory 2 15:30-17:30h
- 31 October 2022 UMinho Ed.1-1.18
- Laboratory 1 16:30-19:00h
- 2 November 2022 INL Lab
- Laboratory 2 16:30-19:00h
- 2 November 2022 INL Lab

Student Feedback




Mathematics for signal processing of scientific instrumentation

| Difficulty level of class? | | | | |
|----------------------------|----------|------------------|---------------|--|
| | Too easy | Just about right | Too difficult | Main comments |
| | 0 | 15 | 5 | More exercise |
| | | | | Unfamiliar with topics, need more introduction |
| | | | | |
| Interesting? | | | | |
| | Yes | Somewhat | No | |
| Too difficult | 2 | 3 | 0 | |
| Just about right | 8 | 7 | 0 | |
| | 10 | 10 | 0 | |
| | | | | |
| Do you want lab class? | | | | |
| | Yes | No | | |
| | 20 | 0 | | |
| | | | | |

Bio-AFM Imaging Modes

INTERNATIONAL IBERIAN NANOTECHNOLOGY LABORATORY

A novel approach to study biological structures




Source: Dufrêne, Y., Ando, T., Garcia, R. et al. Imaging modes of atomic force microscopy for application in molecular and cell biology. *Nature Nanotech* **12**, 295–307 (2017).

Multi-frequency AFM

Cantilever motion






Scheme of the deflection of the cantilever in bimodal AFM. b, Two eigenmodes of the cantilever are excited and detected. Observables associated with both eigenmodes are recorded to determine sample properties such as flexibility, deformation and viscosity.


Multi-frequency AFM

Characterizing ferritin proteins





Ferritin is a universal intracellular <u>protein</u> that stores <u>iron</u> and releases it in a controlled fashion. The protein is produced by almost all living organisms, including archaea, bacteria, algae, higher plants, and animals. It is the primary *intracellular iron-storage protein* in both <u>prokaryotes</u> and <u>eukaryotes</u>, keeping iron in a soluble and non-toxic form. In humans, it acts as a buffer against <u>iron deficiency</u> and <u>iron overload</u>.




Multi-frequency AFM


Characterizing GRoEL proteins




GroEL is a protein which belongs to the <u>chaperonin</u> family of <u>molecular chaperones</u>, and is found in many bacteria. It is required for the proper <u>folding</u> of many proteins

Multifrequency AFM is able to resolve multiple hydration layers of the GrOEL proteins (red dashed line = protein surface)








The first solvation/hydration shell of a <u>sodium</u> <u>ion</u> dissolved in water

High-Speed (HS) AFM

Imaging 1000x faster than before!





Design and Modeling of a High-Speed AFM-Scanner


Georg Schitter, *Member, IEEE*, Karl J. Åström, *Fellow, IEEE*, Barry E. DeMartini, *Member, IEEE*, Philipp J. Thurner, Kimberly L. Turner, and Paul K. Hansma


High-speed AFM

Myosin V walking along actin filaments



Myosins are a <u>superfamily</u> of <u>motor proteins</u> best known for their roles in <u>muscle contraction</u> and in a wide range of other <u>motility</u> processes in <u>eukaryotes</u>. They are <u>ATP</u>-dependent and responsible for <u>actin</u>-based motility






Myosin V walking unidirectionally along an actin filament, showing forward rotation of the leading lever-arm on trailing head detachment from actin

High-speed AFM Imaging example Motorless F1-ATPase undergoing conformational changes



F-ATPase, also known as **F-Type ATPase**, is an <u>ATPase/synthase</u> found in bacterial <u>plasma membranes</u>, in <u>mitochondrial inner membranes</u> (in <u>oxidative phosphorylation</u>, where it is known as Complex V), and in <u>chloroplast thylakoid membranes</u>. It uses a <u>proton</u> gradient to drive ATP synthesis by allowing the passive flux of protons across the membrane down their electrochemical gradient and using the energy released by the transport reaction to release newly formed <u>ATP</u> from the active site of F-ATPase



Red circles indicate the highest positions of the topographs. Since a nucleotide-free β -subunit protrudes higher than ADP- and ATP-bound ones, it is observed that the unbound state rotates anticlockwise.

High-speed AFM Imaging example ESCRT-III protein Snf7 polymerization into spirals on a SLB




The <u>endosomal</u> sorting complexes required for transport (ESCRT) machinery is made up of <u>cytosolic</u> protein complexes, known as ESCRT-0, ESCRT-I, ESCRT-II, and ESCRT-III.

The ESCRT machinery plays a vital role in a number of cellular processes including


- <u>multivesicular body</u> (MVB) biogenesis
- cellular <u>abscission</u>
- viral budding

Snf7 acts a component of the ESCRT-III complex required for the sorting and concentration of proteins resulting in the entry of these proteins into the invaginating vesicles of the multivesicular body (MVB)



High-speed AFM Imaging example MVB biogenesis






PALM/STORM

Maxima localization microscopy





Proposed method for molecular optical imaging


February 1, 1995 / Vol. 20, No. 3 / OPTICS LETTERS 237

E. Betzig



Translation of PALM/STORM concept to the AFM world






Localization atomic force microscopy


<u>George R. Heath, Ekaterina Kots, Janice L. Robertson, Shifra Lansky, George Khelashvili, Harel Weinstein & Simon Scheuring</u>

☐

Data acquisition, analysis and representation











Surface structure, 9 AFM traces & Average-AFM trace















Aquaporin Z (AqpZ) is an integral membrane protein that facilitates transport of water across E. coli cells with a high rate




Annexin A5 (or annexin V) is a cellular protein in the annexin group. Annexin V is commonly used to detect apoptotic cells by its ability to bind to phosphatidylserine, a marker of apoptosis when it is on the outer leaflet of the plasma membrane.




Opto-AFM

Concept & fluorescence laser scanning confocal microscopy / AFM






Laser Scanning Confocal microscopy (LSCM) is most widely applied Opto-AFM system integration:



Jabłoński diagram









Combined imaging of macrophages and fungi





Fluorescence image (left) and correlative AFM images (right) of a macrophage (green) incubated for 3 h with cells from Candida Albicans (blue). AFM images are enlarged views of the dashed areas shown in the fluorescence image. Internalized (bottom) and externalized (top) hyphae featuring major structural differences






CLSM / AFM limitations



Artefacts restricting simultaneous and colocalized data acquistion



| Typical scanning for live cells | Field of View | Image Rep Rate |
|---------------------------------|---------------|----------------|
| CLSM | ~220µm | 1Hz |
| AFM | 100µm | 0.05Hz |


time

- Fluorescence excitation light interacting with the cantilever causes uncontrolled cantilever bending and localized heating due to absorption, mainly of the Au cantilever surface coating
- Mechanical noise transfer from optical microscopy system to AFM camera fans & other moving parts

Simultaneous combined imaging


Differential Spinning Disk (DSD) microscopy








Some key developments in combined AFM and Super-Resolution optical schemes



Stimulated Emission Depletion Microscopy

Concept of STED microscopy






Super-resolution Structured Illumination Microscopy




SR-SIM concept



Fluorescence super-resolution

Comparison of resolution limits







Comparison of spatial resolutions achievable with current commercial super-resolution imaging systems (adapted from Schermelleh et al, 2010).

Super-Resolved Structured Illumination (SR-SIM) and AFM SR-SIM/AFM platform








Photosynthesis adaptation to light environment







Liquid mode AFM hubs Who's who in liquid-mode AFM






Bio-AFM and SARS-CoV-2

Recent work from Alsteens group





Molecular insights into receptor binding energetics and neutralization of SARS-CoV-2 variants

Melanie Koehler, Ankita Ray, Rodrigo A. Moreira, Blinera Juniku, Adolfo B. Poma [™] & David Alsteens [™]

Bio-AFM and AMR

Recent work from Müller/Hiller group






Fig. 7 | **Model of the mechanism of action of polymyxins. 1** Polymyxins initially localize on the OM membrane of Gram-negative bacteria through electrostatic interactions with LPS. **2** On the membrane, polymyxins form hexagonal crystalline structures that increase the membrane surface area, decrease the membrane

bilayer height, and stiffen the membrane. The altered mechanical properties of the OM lead to (3a) the formation of membrane bulging protrusions and (3b) membrane rupture, content leak, which allows progression to the inner membrane and subsequent bacterial death.

How coronavirus enters human cells

INTERNATIONAL IBERIAN NANOTECHNOLOGY LABORATORY

An example case of using ChimeraX to visualize proteins



UCSF ChimeraX

UCSF ChimeraX (or simply ChimeraX) is the next-generation molecular visualization program from the Resource for Biocomputing, Visualization, and Informatics (RBVI), following UCSF Chimera. ChimeraX can be downloaded free of charge for academic, government, nonprofit, and personal use. Commercial users, please see ChimeraX commercial licensing.